BRAUN

No touch + forehead
thermometer

NTF 3000

English 1
Español 17

Contact Us

Questions or Comments? ☎

Call: 1-800-327-7226

eMail: consumerrelations@kaz.com

Visit: www.kaz.com/braun

Write: Consumer Relations, Kaz USA, Inc.,

250 Turnpike Road, Southborough, MA 01772, USA

This product is manufactured by Kaz USA, Inc., under license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

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Visite: www.kaz.com/braun

Escriba a: Consumer Relations, Kaz USA, Inc.,

250 Turnpike Road, Southborough, MA 01772, EEUU

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 Kaz USA, Inc., 250 Turnpike Rd., Southborough, MA 01772, USA

Patents/Patentes: www.kaz.com/patents/braun

Made in China / Hecho en China

© 2014 Kaz USA, Inc. All rights reserved.

Thank you for purchasing the Braun No touch + forehead thermometer (NTF 3000). This thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, the product can provide a stable, accurate reading with each measurement. The thermometer performs a self-test every time it is switched on to always guarantee the accuracy of measurements. The No touch + forehead thermometer is intended for the intermittent monitoring of human body temperature in the home. It is intended for use on people of all ages except pre-term babies or light-for-date babies.

Please read these instructions carefully before using this product and keep the instructions and the thermometer in a safe place.

WARNINGS AND PRECAUTIONS

Keep out of reach of children under 12 years. Never use the thermometer for purposes other than those for which it has been intended. Please follow the general safety precautions when using on children.

Never immerse the thermometer into water or other liquids (not waterproof). Please follow the instructions in the "Care and cleaning" section.

Do not store this thermometer in temperature extremes below -13 °F (-25 °C) or over 140 °F (60 °C) or in excessive humidity (above 95% non-condensing relative humidity). **If thermometer is stored in a location that is cooler or warmer than where it is being used, let it sit in the patient's room for 10 minutes before taking a measurement.**

Do not use the thermometer if there are signs of damage on the scanner or on the thermometer itself. If damaged, do not attempt to repair the product. Never insert a sharp object into the scanner area or any other open surface on the thermometer.

This thermometer consists of high-quality precision parts. Do not drop the device. Protect it from severe impact and shock. Do not twist the device or the sensor.

This thermometer is intended for household use only. Use of this thermometer is not intended as a substitute for consultation with your physician.

Temperature elevation may signal a serious illness, especially in newborns and infants, or in adults who are old, frail, or have a weakened immune system¹. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature on:

- Neonates and infants under 3 months (Consult your physician immediately if the temperature exceeds 99.4°F)
- Patients Over 60 years of age (Fever may be blunted or absent in older patients)
- Patients having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer treatment, chronic steroid treatment, post-surgery)
- Patients who are bedridden (e.g., nursing home patient, stroke, chronic illness, recovering from surgery)
- A transplant patient (e.g., liver, heart, lung, kidney)

This thermometer is not intended for pre-term babies or light-for-date babies.

This thermometer is not intended to interpret hypothermic temperatures. Do not allow children below 12 years to take their temperatures unattended.

Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.


Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on medication should not be assessed solely on temperature readings to determine the severity of their illness.

Do not modify this equipment without the authorization of the manufacturer.

1. See page 15

Product description

- | | |
|----------------------|----------------------|
| 1 Power button | 5 Silent mode switch |
| 2 Temperature button | 6 Battery door |
| 3 LCD display | 7 Guidance light |
| 4 Scanner | 8 Protective cap |

Why Braun No touch + forehead thermometer?

Measurement in under 2 seconds

The innovative infrared technology measures forehead temperature in under 2 seconds, whether you use the touch or no-touch option.


Accurate and reliable

Due to the unique technology, the Braun No touch + forehead thermometer captures the heat naturally given off by the forehead to calculate temperature. Get the same professional accuracy whether touching your child's forehead with the thermometer or holding it up to 2 inches (5 cm) away.

Easy to use

The Braun No touch + forehead thermometer is non-invasive. A measurement can be taken even while a child is sleeping.

The Braun No touch + forehead thermometer is less threatening to a child than a rectal thermometer and easier to use than other methods.

Safe and hygienic

No-touch option helps minimize spreading of germs.

Safe for use on children and adults.

How does Braun No touch + forehead work?

The Braun No touch + forehead thermometer measures infrared energy radiated from the skin at the center of the forehead area. This captured energy (which is twice as much thermal energy compared to a traditional forehead thermometer²), is collected through the sensor and converted to an oral equivalent temperature value.

The Braun No touch + forehead thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

2. A traditional forehead product without any optical system to capture radiated heat.

Temperature guidance feature

Temperature guidance helps you to better understand the meaning of your child's temperature with the color indicated on the display. The screen displays a green, yellow, or red backlight depending on the temperature reading. 10 audible beeps indicate temperatures above 99.4F (37.4C).

It is important to note that a patient's temperature may be impacted by age, gender or medical condition, or by factors such as ambient temperature, time of day, medications taken or the temperature measurement site³.


Color range	Reading
Green	>96.3 - 99.4 °F (>35.7 - 37.4 °C)
Yellow	>99.4 - 101.3 °F (>37.4 - 38.5 °C)
Red	>101.3 - 108.0 °F (>38.5 - 42.2 °C)

NOTE: Backlight remains OFF when below temperature readings are registered: 93.9°F ≤ 96.3°F (34.4 °C ≤ 35.7°C)

3. Herzog L, Phillips SG. Addressing concerns about fever. Clin Pediatr (Phila) 2011 May;50(5):383-90.


How to use your Braun No touch + forehead

1 Remove cap


2 Power on

Press and release the power button once. Backlight will come on and the start up sequence starts.


Ready

When the device is ready and correctly positioned, a horizontal line of dashes ("---") will appear on the screen.

3 Position

Position device on or up to 2 inches (5 cm) away from the center of the forehead, just between the eyebrows. For No touch readings, the yellow guidance light will show you where you are aiming. If the eyebrow area is covered with hair, sweat or dirt, please clean the area and then wait 10 minutes to improve the reading accuracy.

It is important to hold the thermometer and the forehead steady during measurement. Movement will impact the temperature reading.

4 Take temperature

Press the temperature button (you can press and hold the button or press and release it). When the device is placed correctly, the screen displays a dashed line animation while it takes a reading. After the animation sequence (under 2 seconds), the display shows the temperature. The appropriate backlight color is displayed on the screen and the confirmation beep is heard.


If the device is positioned too far away from the forehead, it will prompt you to move closer by displaying the letters "FWD." Slowly move the device toward the forehead until the dashed line animation starts and a reading is displayed.

5 Read temperature

Remove the device and read the temperature.

For no backlight and green backlight readings, you will hear a single long beep for 2 seconds. For yellow and red backlight readings, you will hear 10 short beeps.

To repeat

Go to step 3.

6 To turn off

Press the power button to turn off. Device will also shut off automatically after 60 seconds of no use.

Temperature taking hints

It is important to know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. To determine normal temperature, take multiple readings when healthy. Re-measure with a standard digital thermometer for confirmation. This thermometer displays an oral equivalent temperature reading. Oral readings are approximately 0.9°F (0.5 °C) lower than a rectal digital measurement⁴).

Patient must be inside for 30 minutes before taking a measurement.
Note: Patient and the thermometer should be in the same ambient temperature for at least 10 minutes before taking a reading.

Always hold the thermometer and the forehead steady when taking a reading. Do not move the thermometer until you hear the final beep.

For patients measuring their own temperature, use the "touch" option instead of "no touch."

Don't take a measurement while or immediately after nursing a baby.

Patients should not drink, eat, or be physically active before/while taking the measurement. Remove hats and wait 10 minutes before taking a reading.

Remove dirt or hair from the forehead before taking a measurement. Bangs may cause higher readings. Wait at least 10 minutes after washing the forehead area before taking a reading.

Always take the temperature exactly as directed. Temperature results may vary if positioned in the wrong location.


In the following situations it is recommended that three temperatures in the same location be taken and the highest one taken as the reading:

- Newborn infants in the first 100 days.
- Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
- When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.


Changing the temperature scale

- 1 Open battery compartment and remove the batteries.

- 2 °C / °F switch is found in the center of the battery compartment.

- 3 Slide switch to °C or °F to set preferred temperature scale.


- 4 Replace the batteries and close the battery door.

4. Sund-Levander M, Forsberg C, Wahren LK. Normal oral, rectal, tympanic and axillary body temperature in adult men and women: a systematic literature review. Scand J Caring Sci 2002 June;16(2):122-8
Herzog L, Phillips SG. Addressing concerns about fever. Clin Pediatr (Phila) 2011 May;50(5):383-90.

Changing the sound mode

The Braun No touch + forehead thermometer allows you to silence the beeps on the thermometer to avoid waking a sleeping child. To use this feature, slide the switch to the silent mode position and take a temperature reading. An icon will appear on the screen to show that the thermometer is in silent mode.


Calibration

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instructions, periodic re-adjustment is not required. If at any time you question the accuracy, please contact Consumer Relations.

Manufacturing date is given by the LOT number located inside the battery compartment. The first three (3) digits after LOT represent the Julian date that the product was manufactured and the next two (2) digits represent the last two numbers of the calendar year the product was manufactured. The last identifiers are the letters that represent the manufacturer.

An example: LOT 11614tav, this product has been manufactured on the 11th day of the year 2014.

Errors and troubleshooting

Error Message	Situation	Solution
	When 20% of battery life is left, the display flashes the low battery warning symbol. The device can still work until the battery life has 0% left.	Replace battery.
	If the steady battery icon is the only symbol shown on the display, the device cannot work. The battery should be replaced immediately.	Replace battery.
	This message displays when measured temperature is lower than 93.9 °F (34.4 °C) or higher than 108 °F (42.2 °C) or when the room temperature is outside the operating range of 59 °F-104 °F (15°C-40°C).	Make sure the protective cap is removed. Re-measure the temperature again, carefully following the instruction in "How to use" section.
	Blank display. Thermometer does not have power.	Please check if the battery has been loaded correctly. Also check polarity (<+> and <->) of batteries. Contact Consumer Relations if thermometer still does not function.

Care and cleaning

The thermometer should be cleaned in between uses, as necessary. Use an alcohol swab or cotton swab moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never soak the device in water or other cleaning liquids. Wait 10 minutes after cleaning before taking a temperature measurement. Make sure to replace the protective cap after use to prevent scratches or damage to the thermometer.

Never insert a sharp object into the scanner area or any other open surface on the thermometer.

Clean sensor mirror before each use or as needed.

Replacing the batteries

The Braun No touch + forehead thermometer comes with 2 AA batteries. Replace with 2 new AA batteries when the flashing battery symbol appears on the LCD display. To change the batteries, slide open the battery cover and remove batteries. Replace the batteries being sure to align properly as indicated inside the battery compartment. Remove the battery from the product if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.

To protect the environment, dispose of the product and empty batteries at appropriate collection sites according to national or local regulations.

Limited three year warranty

Kaz warrants this product against any defects that are due to faulty material or workmanship for a period of three years from the original date of consumer purchase or receipt as a gift. This warranty applies when used for normal household use in accordance with the operating instruction manual and excludes the battery and damage to the product resulting from accident or misuse. This product is not warranted when used in a professional environment.

In no event shall Kaz be liable for any special, incidental, indirect or consequential damages in connection with the purchase or use of this product or costs over the original cost of the product.

If the product should not perform to specifications within the warranty period, please contact Consumer Relations. Kaz reserves the right to replace a non conforming product with the most comparable product currently available.

While this warranty gives you specific legal rights, you may also have other rights which vary from state to state.

Product specification

Type:	No touch + forehead thermometer (NTF 3000)
Measuring range:	93.9 °F (34.4 °C) – 108 °F (42.2 °C)
Resolution:	0.1 °F (0.1 °C)
Laboratory Accuracy	± 0.4°F for the range 95°-107.6°F (± 0.2°C for 35-42°C) (Ambient Temperature: 59°F to 104°F (15°C to 40°C)
Display:	Liquid Crystal Display, 4 digits plus special icons
Acoustic:	Green backlight = single long beep for 2 sec duration
Audio:	Red or Yellow backlight: 10 short beeps
Operating temperature:	59 °F – 104 °F (15 °C – 40 °C)
Automatic Switch-off:	Approx. 60 seconds after last measurement has been taken
Weight:	3.51oz. (with battery), 2.72oz. (w/o battery)

Long term storage ranges

Storage/transport temperature: -13 °F to 140 °F (-25 °C to 60 °C)
 Humidity: 15-95% non condensing
 Battery: (2) AA Batteries - at least 1000 measurements
 Pressure: 700-1060hPA (0.7-1.06 atm)

This infrared thermometer meets requirements established in ASTM Standard E 1965-98 for the thermometer system. Full responsibility for the conformance of the product to the standard is assumed by Kaz USA, Inc., 250 Turnpike Rd., Southborough, MA 01772, USA.

ASTM laboratory accuracy requirements for the thermometer only in the display range of 96.8 °F to 102.2 °F (36 °C to 39 °C) for infrared thermometers is ± 0.4 °F (± 0.2 °C), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is ± 0.2 °F (± 0.1 °C).

This device conforms to the following standards:

IEC 60601-1:2006 Medical electrical equipment. General requirements for basic safety and essential performance

ISO 80601-2-56 Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement

ASTM 1965E – ASTM Standard for Infrared Clinical Thermometers

EN 12470-5: 2003 «Clinical thermometers» – Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.

IEC 60601-1-2: General requirements for basic safety and essential performance


Collateral Standard: Electromagnetic compatibility – Requirement and tests

NOTE: Do not use this device in the presence of electromagnetic or other interference outside the normal range specified in IEC 60601-1-2.

EN 980: 2008 Symbols for use in labeling of medical devices.

EN 1041: 2008 Information supplied by the manufacturer of medical devices.

EN 60601-1-11: 2010 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

**TYPE BF EQUIPMENT**

Internally Powered ME Equipment

Continuous Operation

IP20: Protected against solid foreign objects of 12.5 mm diameter and greater.

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC. For detailed description of EMC requirements please contact Consumer Relations.

Portable and mobile RF communications equipment can affect **MEDICAL ELECTRICAL EQUIPMENT**.

- Dinarello CA, Porat R. Fever and hyperthermia. In: Fauci AS, Braunwald E, Hauser SL, Longo DL, Loscalzo J, editors. Harrison's Internal Medicine. 17th ed. New York, NY: McGraw Hill Medical; 2008.
- American College of Emergency Physicians. Clinical policy for children younger than 3 years presenting to the emergency department with fever. Ann Emerg Med 2003;43(4):530-43.
- Byington CL, Rittichier KK, Bassett KE et al. Serious bacterial infections in febrile infants younger than 90 days of age: the importance of ampicillin-resistant pathogens. Pediatrics 2003 May;111(5 Pt 1):964-8.
- High KP, Bradley SF, Gravenstein S et al. Clinical practice guideline for the evaluation of fever and infection in older adult residents of long-term care facilities: 2008 update by the Infectious Diseases Society of America. J Am Geriatr Soc 2009 March;57(3):375-94.
- Fox RH, Woodward PM, Exton-Smith AN, Green MF, Donnison DV, Wicks MH. Body temperatures in the elderly: a national study of physiological, social, and environmental conditions. Br Med J 1973 January 27;1(5847):200-6.
- Holt RIG, Cockram C, Flyvbjerg A, Goldstein BJ. Textbook of Diabetes. Wiley-Blackwell; 2010
- Doan T, Melvold R., Viselli S., Waltenbaugh C. Lippincott's Illustrated Reviews: Immunology, Second Edition. Lippincott Williams & Wilkins; 2013

Non-Life Support Equipment Separation Distance Calculation (3Vrms / 3V/m compliance)			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz in ISM bands $d = \sqrt{\frac{3.5}{V_1}} \sqrt{P}$	80 MHz to 800 MHz $d = \sqrt{\frac{3.5}{E_1}} \sqrt{P}$	800 MHz to 2.5 GHz $d = \sqrt{\frac{7}{E_1}} \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

Guidance and manufacturer's declaration – electromagnetic immunity			
The NTF 3000 is intended for use in the electromagnetic environment specified below. The customer or the user of the NTF 3000 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV Contact ±8kV Air	Complies	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5GHz	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	Not Applicable (no electrical cabling)	Separation distance calculation provided below. If a known transmitter is present the specific distance can be calculated using the equations.
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	Not Applicable	The ME equipment is solely battery powered.
Surge IEC 61000-4-5	±1kV differential ±2kV common	Not Applicable	
Power frequency magnetic field IEC 61000-4-8	3 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% dip 0.5 cycle 60% dip 5 cycles 70% dip 25 cycles 95% dip 5 sec.	Not Applicable	The ME equipment is solely battery powered.

Guidance and manufacturer's declaration – electromagnetic emissions		
The NTF 3000 equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the NTF 3000 should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	Complies
Harmonic emissions IEC 61000-3-2	Not Applicable	The ME equipment is solely battery powered.
Voltage fluctuations/flicker emissions	Not Applicable	

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