Go Relax™ Lumbar Cushion

Cojín de Masaje
Coussin de Massage

FEATURES

- Six accupressure nodes
- Pressure activated switch
- Built-in high density foam cushion

OPERATION

1. Slide cover off battery compartment, located beneath zipper on backside of cushion. Install 2 AA batteries.
2. Lean against cushion to activate vibrating massage.
3. Unit will shut off automatically upon release of pressure.

CARE AND MAINTENANCE

- Wipe clean with a soft, damp cloth.
- Never submerge the unit in any liquid.
- Store in a dry and cool place. Avoid contact with sharp edges or pointed objects which might cut or puncture the surface.

NOTE

Accupressure therapy is not meant as a substitute for any medical treatment prescribed by doctors and health care providers, but may be used in addition to other treatments. Wagan Corp. makes no claim to diagnose, prevent, treat or cure disease.

WAGAN Corp. Limited Warranty

The WAGAN Corp. warranty is limited to products sold only in the United States.

All Wagan Tech products must be registered within (30) days of purchase to activate its warranty. To register your product, please visit <http://tinyurl.com/waganwarranty>. Be sure to keep the original receipt as it will be required when returning a product under the warranty.

Warranty Duration: This product is warranted to the original purchaser for a period of one (1) year from the original purchase date, to be free of defects in material and workmanship. WAGAN Corporation disclaims any liability for consequential damages. In no event will WAGAN Corporation be responsible for any amount of damages beyond the amount paid for the product at retail.

Warranty Performance: During the above one (1) year warranty period, a product with a defect will be replaced with a comparable model when the product is returned to WAGAN Corporation with an original store receipt. The replacement product will be in warranty for the balance of the original one (1) year warranty period.

To return a defective item, please contact WAGAN Corporation at (800) 231-5806 to obtain a Returned Merchandise Authorization number (RMA#), and return instructions. Each item returned will require a separate RMA#. After you have received the RMA# and the return instructions from WAGAN Corporation, please follow the instructions and send the item with PREPAID SHIPPING, along with all of the required documentation, a complete explanation of the problem, your name, address and daytime phone number. WAGAN Corporation will, at its option, replace or repair the defective part.

A Returned Merchandise Authorization number (RMA#) is REQUIRED when sending in any defective item. WAGAN Corporation is not responsible for any item(s) returned without an official Returned Merchandise Authorization number. The item(s) must be returned with prepaid shipping. WAGAN Corporation is not responsible for any shipping charges incurred in returning the item(s) back to the company for repair or replacement. This warranty is void if the product has been damaged by accident, in shipment, unreasonable use, misuse, neglect, improper service, commercial use, repairs by unauthorized personnel or other causes not arising out of defects in materials or workmanship. This warranty does not extend to any units which have been used in violation of written instructions furnished.

Warranty Disclaimers: This warranty is in lieu of all warranties expressed or implied and no representative or person is authorized to assume any other liability in connection with the sale of our products. There shall be no claims for defects or failure of performance or product failure under any theory of tort, contract or commercial law including, but not limited to negligence, gross negligence, strict liability, breach of warranty and breach of contract.

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©2014 Wagan Corporation. All Rights Reserved
Wagan Tech and *wagan.com* are trademarks used by Wagan Corporation
Garantie Limitée de WAGAN Corporation
La Garantie WAGAN Corporation est limitée aux produits vendus sur le territoire des États-Unis.
Hayward, CA 94544, U.S.A.
© 2014 Wagan Corporation. Tous droits réservés.